



PRIMA DENTAL


PREDATOR®

Metal Cutting Carbide Burs

Predator® Turbo and Predator® Accu-Prep burs cut faster than straight bladed burs and smoother than regular dentated burs.

They are the strongest, sharpest and most durable burs available with unparalleled concentricity. They are manufactured using the latest state of the art Swiss-made precision machine tools.

This technology lifts the Predator® burs to the height of premium dental instruments and dentistry.



METAL CUTTING BURS



PREDATOR®
TURBO
PR 3T

US NO.	PR 3T
DIAMETER 1/10 MM	012
H/LENGTH (MM)	4.0

Single piece carbides



PREDATOR®
DX

US NO.	PR DX
DIAMETER 1/10 MM	012
H/LENGTH (MM)	4.0

For separating
porcelain fused
to single piece
metal crowns



PREDATOR®
TURBO

US NO.	PR 2T	PR 1557T	PR 1T	PR 330T	PR 245T	PR 557T	PR 558T	PR35	PR 701T	PR 702T
FG	012	010	010	008	009	010	012	010	12.0	16.0
FG SURG	012		010			010				
H/LENGTH (MM)	4.0	3.9	1.9	2.0	2.7	3.7	3.8	0.79	3.80	4.30

Single piece carbides



PREDATOR®
TURBO

US NO.	PR 2R	PR 4R	PR 6R	PR 8R
FG	010	014	018	023
FG SURG	010	014	018	023
H/LENGTH (MM)	0.83	1.23	1.57	2.05



PREDATOR® S1

US NO.	ISO NO.	S1
FG (19MM)	500 314 137 006	012

Single piece carbide, domed
Fissure Plain Cut



PREDATOR® ACCU-PREP

US NO.	379	847	856	856	856
DIAMETER 1/10 MM	023	016	010	016	018
H/LENGTH (MM)	4.3	8.0	8.0	8.0	8.0

Carbides for crown preparation



Shoulder Preps



Chamfer Preps



Lingual / Occlusal



To explore
our Predator®
Diamond range,
visit our website
or get in touch

www.primadental.com

PRIMA DENTAL

Stephenson Drive, Gloucester,
GL2 2HA, United Kingdom

T +44 1452 729751

F +44 1452 307187

E sales@primadental.com

W www.primadental.com



PRIMA DENTAL